

Product Datasheet

Anti-h PIVKA-II 12102 SPTN-5 100871

Product Name	Anti-h PIVKA-II 12102 SPTN-5
Catalog Number	100871
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA, CLIA
Alternative Names	Protein Induced by Vitamin K absence or Antagonist-II
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Amino acid region 11-28 (GNLERECVEETCSYEEAF)
Purity	≥ 95 %
Affinity constant	N/D
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-7.1
Cross Reactivity	Does not recognize Prothrombin.
Specificity	Antibody recognizes human PIVKA-II (Protein Induced Vitamin K absence or Antagonist-II)
Matched Pairs	PIVKA-II 12103 and PIVKA-II 12106
Shelf Life	Unspecified