

Product Datasheet

Anti-h PIVKA-II 12102 SPTN-5 100871

| | |
|----------------------------|--|
| Product Name | Anti-h PIVKA-II 12102 SPTN-5 |
| Catalog Number | 100871 |
| Description | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Species/Host | Mouse |
| Tested Applications | CLIA, FIA |
| Alternative Names | Protein Induced by Vitamin K absence or Antagonist-II |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/- 10 %) |
| Storage | +2-8 °C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Epitope | Amino acid region 11-28 (GNLERECVEETCSYEEAF) |
| Purity | ≥ 95 % |
| Affinity constant | N/D |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| IEF Profile | 6.3-7.1 |
| Cross Reactivity | Does not recognize Prothrombin. |
| Specificity | Antibody recognizes human PIVKA-II (Protein Induced Vitamin K absence or Antagonist-II) |
| Matched Pairs | PIVKA-II 12103 and PIVKA-II 12106 |
| Shelf Life | Unspecified |