

Product Datasheet

Monkeypox Virus (A5L & A11L) Antigen, Recombinant >85% LA804

Product Name Monkeypox Virus (A5L & A11L) Antigen, Recombinant >85%

Catalog Number LA804

Brand EastCoast Bio

Form/Appearance Affinity Purified Recombinant Antigen (His-tag)

Storage ≤ -20°C. Avoid repeated freezing and thawing.

Epitope A5L & A11L proteins

Purity >85%

Buffer 50 mM Tris-HCl, pH 8.0 with 100 mM NaCl

Expression System E. coli

Preservatives None