

Product Datasheet

Monkeypox Virus (H3L) Antigen, Recombinant 85% LA805

Product Name	Monkeypox Virus (H3L) Antigen, Recombinant 85%
Catalog Number	LA805
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	≤ -20°C. Avoid repeated freezing and thawing.
Epitope	H3L protein
Purity	>85%
Buffer	50 mM Tris-HCl, pH 8.0 with 100 mM NaCl
Expression System	E. coli
Preservatives	None