

Product Datasheet

Dinitrophenyl (DNP) Antibody HM1328

Product Name	Dinitrophenyl (DNP) Antibody
Catalog Number	HM1328
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	DNP-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	PBS, pH 7.0
Preservatives	0.02% Sodium Azide