

## Product Datasheet

### Dinitrophenyl (DNP) Antibody HM1328

<b>Product Name</b>	Dinitrophenyl (DNP) Antibody
<b>Catalog Number</b>	HM1328
<b>Description</b>	Ion Exchange Purified Monoclonal Antibody
<b>Immunogen</b>	DNP-KLH
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified monoclonal antibody
<b>Storage</b>	Short term: 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	Mouse IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 % by SDS-PAGE
<b>Buffer</b>	Phosphate-buffered saline, pH 7.0 ± 0.2, 0.02% sodium azide as a preservative.
<b>Specificity</b>	Dinitrophenyl (DNP)