

Product Datasheet

Methaqualone-BSA Antigen 170072

Product Name	Methaqualone-BSA Antigen
Catalog Number	170072
Description	Synthetic Methaqualone, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Concentration	Lot-specific
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P82-99-81A. It has been produced using the same conjugation protocol.
Buffer	Distilled Water
Matched Pairs	P02-99-81R-IF, P02-99-82R-IF or P02-99-83R-IF (detection)
Shelf Life	60 months