

# Product Datasheet

## Hydrocodone-BSA Antigen 170054

Product Name	Hydrocodone-BSA Antigen
Catalog Number	170054
Description	Synthetic Hydrocodone, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Concentration	Lot-specific
Storage	Short term storage ( $\leq 2$ weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P83-99-71A and LA289. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	JP101 or JP102 (detection)
Shelf Life	60 months