

Product Datasheet

K2-BSA Antigen 170082

Product Name	K2-BSA Antigen
Catalog Number	170082
Description	Synthetic K2, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Concentration	Lot-specific
Storage	Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P83-99-51A and LA288. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	JP098, JP099 or JP100 (detection)
Shelf Life	60 months