

Product Datasheet

Ethyl glucuronide (EtG) Antibody HM1400

| | |
|----------------------------|--|
| Product Name | Ethyl glucuronide (EtG) Antibody |
| Catalog Number | HM1400 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Goat |
| Tested Applications | LF, ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Clonality | Monoclonal |
| Purity | ≥ 95 % by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 ± 0.2 |
| Specificity | Ethyl glucuronide |
| Matched Pairs | LA809 (capture) |