

## Product Datasheet

### Ethyl glucuronide (EtG) Antibody HM1400

<b>Product Name</b>	Ethyl glucuronide (EtG) Antibody
<b>Catalog Number</b>	HM1400
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 % by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4 ± 0.2
<b>Specificity</b>	Ethyl glucuronide
<b>Matched Pairs</b>	LA809 (capture)