

Product Datasheet

Ethyl glucuronide (EtG)-BSA Antigen LA809L

Product Name	Ethyl glucuronide (EtG)-BSA Antigen
Catalog Number	LA809L
Tested Applications	ELISA, Lateral Flow
Brand	EastCoast Bio
Form/Appearance	Lyophilized; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C.
Note	Avoid repeated freezing and thawing. Pre-lyophilization buffer 10 mM PBS, pH 7.4
Purity	$\geq 90\%$
Matched Pairs	HM1400 (detection)