

## Product Datasheet

### Ethyl glucuronide (EtG)-BSA Antigen LA809L

<b>Product Name</b>	Ethyl glucuronide (EtG)-BSA Antigen
<b>Catalog Number</b>	LA809L
<b>Tested Applications</b>	ELISA, Lateral Flow
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized; Antigen-BSA conjugate
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C.
<b>Note</b>	Avoid repeated freezing and thawing. Pre-lyophilization buffer 10 mM PBS, pH 7.4
<b>Purity</b>	$\geq 90\%$
<b>Matched Pairs</b>	HM1400 (detection)