

## Product Datasheet

### Monkeypox Virus (MPXV-A29L) Antibody HM1338

<b>Product Name</b>	Monkeypox Virus (MPXV-A29L) Antibody
<b>Catalog Number</b>	HM1338
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	in vitro cell culture
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Purified A29L protein
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 % by SDS-PAGE
<b>Buffer</b>	Tris Glycine, pH 7.0-7.4
<b>Matched Pairs</b>	HM1339 (detection)