

Product Datasheet

Monkeypox Virus (MPXV-A29L) Antibody HM1339

Product Name	Monkeypox Virus (MPXV-A29L) Antibody
Catalog Number	HM1339
Description	Protein A Purified Monoclonal Antibody
Species/Host	in vitro cell culture
Tested Applications	ELISA, LF
Immunogen	Purified A29L protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	Tris Glycine, pH 7.0-7.4
Matched Pairs	HM1338 (capture)