

Product Datasheet

Anti-Human CD31 Monoclonal Antibody, Azide Free Clone B-B38 852.560.005

Product Name	Anti-Human CD31 Monoclonal Antibody, Azide Free Clone B-B38
Catalog Number	852.560.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC
Immunogen	Recombinant human CD31
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Platelet Endothelial Cell Adhesion Molecule-1 (PECAM-1), GPIIa, a 130-140 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells