

Product Datasheet

Ubiquitin carboxy terminal Hydrolase (UCHL1) Antibody HM1420

Product Name	Ubiquitin carboxy terminal Hydrolase (UCHL1) Antibody
Catalog Number	HM1420
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Human UCLH1 (HEK293 expressed)
Alternative Names	UCH-L1, PGP9.5, ubiquitin C-terminal hydrolase L1, Traumatic brain injury (TBI)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1415, HM1416, HM1417, HM1418, HM1419 (capture)
Preservatives	0.05% Sodium Azide