

## Product Datasheet

### Ubiquitin carboxy terminal Hydrolase (UCHL1) Antibody HM1422

<b>Product Name</b>	Ubiquitin carboxy terminal Hydrolase (UCHL1) Antibody
<b>Catalog Number</b>	HM1422
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Human UCLH1 (HEK293 expressed)
<b>Alternative Names</b>	UCH-L1, PGP9.5, ubiquitin C-terminal hydrolase L1, Traumatic brain injury (TBI)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 %
<b>Buffer</b>	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
<b>Matched Pairs</b>	HM1415, HM1416, HM1417, HM1418, HM1419 (capture)
<b>Preservatives</b>	0.05% Sodium Azide