

Product Datasheet

Anti-Human CD31 Monoclonal Antibody, FITC Conjugated Clone B-B38 852.561.010

Product Name	Anti-Human CD31 Monoclonal Antibody, FITC Conjugated Clone B-B38
Catalog Number	852.561.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	IHC, FC
Immunogen	Recombinant human CD31
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5 % BSA and 0.1 % sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Platelet Endothelial Cell Adhesion Molecule-1 (PECAM-1), GPIIa, a 130-140 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells