

Product Datasheet

Anti-Human CD48 Monoclonal Antibody, Unconjugated Clone MEM-102 852.613.020

Product Name	Anti-Human CD48 Monoclonal Antibody, Unconjugated Clone MEM-102
Catalog Number	852.613.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IP, WB, IHC, FC
Immunogen	Burkitt's lymphoma cell line Raji
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD48 antigen (Blast-1), a 40-47 kDa GPI-anchored membrane protein (immunoglobulin supergene family) widely expressed on hematopoietic cells