

Product Datasheet

Anti-Human CD50 Monoclonal Antibody,Azide Free Clone B-R1 852.670.005

Product Name Anti-Human CD50 Monoclonal Antibody,Azide Free Clone B-R1

Catalog Number 852.670.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, IP, FC

Immunogen U266 cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Intercellular Adhesion Molecule-3 (ICAM-3), a 110-140 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells