

Product Datasheet

Anti-Human CD54 Monoclonal Antibody, Azide Free Clone B-H17 852.690

Product Name	Anti-Human CD54 Monoclonal Antibody, Azide Free Clone B-H17
Catalog Number	852.690
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC, Functional assay, ELISA, IP
Immunogen	Activated human umbilical vein endothelial (HUVEC)
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Intercellular Adhesion Molecule-1 (ICAM-1), a 75-115 kDa protein
Biological Activity	Inhibits PHA stimulated T-cell proliferation
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells