

Product Datasheet

Anti-Human CD54 Monoclonal Antibody, Azide Free Clone B-H17 852.690.005

Product Name Anti-Human CD54 Monoclonal Antibody, Azide Free Clone B-H17

Catalog Number 852.690.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, IP, Functional assay, FC

Immunogen Activated human umbilical vein endothelial (HUVEC)

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Intercellular Adhesion Molecule-1 (ICAM-1), a 75-115 kDa protein

Biological Activity Inhibits PHA stimulated T-cell proliferation

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells