

Product Datasheet

Anti-Human CD102 Monoclonal Antibody, Azide Free Clone B-T1 852.760.005

Product Name	Anti-Human CD102 Monoclonal Antibody, Azide Free Clone B-T1
Catalog Number	852.760.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, IP, Functional assay, FC
Immunogen	ICAM-2 transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Intercellular Adhesion Molecule-2 (ICAM-2), a 55-65 kDa protein
Biological Activity	Inhibits cell adhesion and T-cell activation
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells