

Product Datasheet

Anti-Human CD106 Monoclonal Antibody, Azide Free Clone B-K9 852.770.005

Product Name	Anti-Human CD106 Monoclonal Antibody, Azide Free Clone B-K9
Catalog Number	852.770.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Functional assay
Immunogen	Activated human umbilical vein endothelial cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Vascular Cellular Adhesion Molecule-1 (VCAM-1), a 110 kDa protein
Biological Activity	Inhibits cell attachment to VCAM-1
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells