

Product Datasheet

Anti-Human CD62L Monoclonal Antibody, PE Conjugated Clone B-S13 852.782.010

Product Name	Anti-Human CD62L Monoclonal Antibody, PE Conjugated Clone B-S13
Catalog Number	852.782.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	CD62L transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Leucocyte Endothelial Cell Adhesion Molecule-1 (LECAM-1), L-selectin, a 74-95 kDa protein. Cross reacts with mouse CD62L protein in cytometry.
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells