

## Product Datasheet

### Anti-h D-Dimer 1402 SPTNZ-5 100205

<b>Product Name</b>	Anti-h D-Dimer 1402 SPTNZ-5
<b>Catalog Number</b>	100205
<b>Description</b>	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
<b>Tested Applications</b>	FIA, IT
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/- 10 %)
<b>Storage</b>	+2-8°C
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Epitope</b>	N/D
<b>Purity</b>	≥ 95 %
<b>Affinity constant</b>	N/A
<b>Associated Products</b>	Native D-Dimer, Lee Biosolutions, Cat. 200-09, 200-12 and 200-13
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.05 % Sulfobetaine, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	6.5-8.4
<b>Cross Reactivity</b>	Recognizes human fibrinogen
<b>Specificity</b>	Antibody recognizes human D-dimer
<b>Shelf Life</b>	12 months