

Product Datasheet

Anti-Human CD107a Monoclonal Antibody, Azide Free Clone B-T47 852.900.005

Product Name	Anti-Human CD107a Monoclonal Antibody, Azide Free Clone B-T47
Catalog Number	852.900.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	U937 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises a 100-120 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells