

# Product Datasheet

## Anti-Human Galectin-10 Monoclonal Antibody, Unconjugated Clone B-F42 852.963.020

Product Name	Anti-Human Galectin-10 Monoclonal Antibody, Unconjugated Clone B-F42
Catalog Number	852.963.020
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IHC, FC
Immunogen	Recombinant human Galectin-10
Alternative Names	Charcot-Leyden Crystal protein   CLC
Brand	Diaclone
Form/Appearance	Liquid
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Note	Galectin-10 is an intracellular antigen. Fixation with paraformaldehyde 4% and permeabilisation with saponin 0.1% are required before staining. Use 10 µl to label 10 <sup>6</sup> cells
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Specificity	Recognises both natural and recombinant human Galectin-10
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells