

## Product Datasheet

### Anti-Human Galectin-10 Monoclonal Antibody, Unconjugated Clone B-F42 852.963.020

<b>Product Name</b>	Anti-Human Galectin-10 Monoclonal Antibody, Unconjugated Clone B-F42
<b>Catalog Number</b>	852.963.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC, IHC
<b>Immunogen</b>	Recombinant human Galectin-10
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
<b>Note</b>	Galectin-10 is an intracellular antigen. Fixation with paraformaldehyde 4% and permeabilisation with saponin 0.1% are required before staining. Use 10 µl to label 10 <sup>6</sup> cells
<b>Isotype</b>	Mouse IgG1 Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognises both natural and recombinant human Galectin-10
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c spleen cells