

Product Datasheet

Anti-Human CD95 Monoclonal Antibody, Azide Free Clone B-G27 853.000.005

Product Name	Anti-Human CD95 Monoclonal Antibody, Azide Free Clone B-G27
Catalog Number	853.000.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Functional assay, IHC, FC
Immunogen	Recombinant Fas antigen
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Fas, APO1 antigen, a 45 kDa protein
Biological Activity	Triggers Fas induced apoptosis
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells