

Product Datasheet

Anti-Human CD178 Monoclonal Antibody, Azide Free Clone B-R17 853.070.005

Product Name Anti-Human CD178 Monoclonal Antibody, Azide Free Clone B-R17

Catalog Number 853.070.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, Functional assay, FC

Immunogen FasL transfected 293T cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises APO1L, Fas ligand (FAS L), a 40 kDa protein

Biological Activity Inhibits soluble CD178 induced apoptosis

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells