

Product Datasheet

Anti-Human CD178 Monoclonal Antibody, Biotin Conjugated Clone B-R17 853.074.010

Product Name	Anti-Human CD178 Monoclonal Antibody, Biotin Conjugated Clone B-R17
Catalog Number	853.074.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	FC, ELISA
Immunogen	FasL transfected 293T cell line
Brand	Diaclone
Form/Appearance	Lyophilised in phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. For longer storage freeze aliquots.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises APO1L, Fas ligand (FAS L), a 40 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells