

Product Datasheet

Anti-Human CD253 Monoclonal Antibody, Biotin Conjugated Clone B-S23 853.084.010

Product Name	Anti-Human CD253 Monoclonal Antibody, Biotin Conjugated Clone B-S23
Catalog Number	853.084.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	Recombinant human TRAIL
Brand	Diaclone
Form/Appearance	Lyophilised in phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. For longer storage freeze aliquots.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the TRAIL, APO2L antigen, a 32 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells