

## Product Datasheet

### Albumin, Human Serum - 25% Solution 101-15

<b>Product Name</b>	Albumin, Human Serum - 25% Solution
<b>Catalog Number</b>	101-15
<b>Alternative Names</b>	Serum ALB, HSA, ALB
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Solution in saline with sodium caprylate and N-acetyl tryptophan as stabilizers.
<b>Concentration</b>	23.5 - 26.5 g/dL (Biuret)
<b>Storage</b>	+2-8 °C
<b>Purity</b>	≥ 96% (Electrophoresis)
<b>CAS Number</b>	70024-90-7
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	66,500
<b>Shelf Life</b>	5 years