

## Product Datasheet

### Aprotinin 125-09

<b>Product Name</b>	Aprotinin
<b>Catalog Number</b>	125-09
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Clear, colorless solution
<b>Concentration</b>	1-3 mg protein/mL (Biuret)
<b>Storage</b>	+2-8 °C
<b>CAS Number</b>	9087-70-1
<b>Source</b>	Bovine Lung
<b>Molecular Weight</b>	6511.44
<b>Shelf Life</b>	2 years
<b>Biological Activity</b>	6 - 16 TIU/mL
<b>Specific Activity</b>	≥ 5 TIU/mg protein