

Product Datasheet

Aprotinin 125-09

Product Name	Aprotinin
Catalog Number	125-09
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear, colorless solution
Concentration	1-3 mg protein/mL (Biuret)
Storage	2-8° C
CAS Number	9087-70-1
Source	Bovine Lung
Molecular Weight	6511.44
Shelf Life	2 years
Biological Activity	6 - 16 TIU/mL
Specific Activity	≥ 5 TIU/mg protein