

Product Datasheet

Apolipoprotein CIII (ApoCIII), Human - Ultra Pure 125-21

Product Name	Apolipoprotein CIII (ApoCIII), Human - Ultra Pure
Catalog Number	125-21
Alternative Names	ApoCIII, ApoC3
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	Typically 1.0 mg/mL (Lowry)
Storage	-20 °C
Purity	≥ 98% (SDS-PAGE)
Source	Human Plasma
Shelf Life	2 years