

Product Datasheet

Apolipoprotein E (ApoE), Human - High Purity 125-29

Product Name	Apolipoprotein E (ApoE), Human - High Purity
Catalog Number	125-29
Alternative Names	ApoE, Apo E
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	0.5 - 5.0 mg/mL ($E^{0.1\%}_{280nm} = 1.31$)
Storage	-60 °C or below
Purity	≥ 95% (SDS-PAGE)
Source	Human Plasma
Molecular Weight	34000
Shelf Life	2 years
Expression System	E. coli