

## Product Datasheet

### Anti-Human CD1a Monoclonal Antibody, Unconjugated Clone HI149 853.153.020

<b>Product Name</b>	Anti-Human CD1a Monoclonal Antibody, Unconjugated Clone HI149
<b>Catalog Number</b>	853.153.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC, IHC
<b>Immunogen</b>	Human thymocytes
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Concentration of 1mg/1ml in phosphate buffered saline containing 15mM sodium azide
<b>Storage</b>	Store at +2-8°C. DO NOT FREEZE.
<b>Note</b>	We recommend that each laboratory determine an optimum working titre for use in its particular application
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Reacts with CD1a (T6), a 49 KDa polypeptide associated with b2-microglobulin expressed on cortical thymocytes (strongly), Langerhans cells and dendritic cells