

## Product Datasheet

### Carcinoembryonic Antigen (CEA), Goat Polyclonal anti-Human Antibody 151-12

<b>Product Name</b>	Carcinoembryonic Antigen (CEA), Goat Polyclonal anti-Human Antibody
<b>Catalog Number</b>	151-12
<b>Species/Host</b>	Goat
<b>Immunogen</b>	Purified, native human CEA
<b>Alternative Names</b>	anti-CEA
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	≥ 1 mg/mL
<b>Storage</b>	2-8° C
<b>Purity</b>	High Purity
<b>Source</b>	Polyclonal Antibody
<b>Shelf Life</b>	2 years