

Product Datasheet

Anti-Human CD11b Monoclonal Antibody, FITC Conjugated Clone MEM-174 853.161.010

Product Name	Anti-Human CD11b Monoclonal Antibody, FITC Conjugated Clone MEM-174
Catalog Number	853.161.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Human granulocytes
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) containing 15 mM sodium azide and 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA)
Storage	Store at +2-8°C protected from light. DO NOT FREEZE.
Note	For flow cytometry we recommend using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension
Isotype	IgG2a
Clonality	Monoclonal
Specificity	Recognizes CD11b antigen (MAC-1), a transmembrane protein type 1 mainly expressed on monocytes, granulocytes and NK-cells