

## Product Datasheet

### Anti-Human CD11b Monoclonal Antibody, Unconjugated Clone MEM-174 853.163.020

<b>Product Name</b>	Anti-Human CD11b Monoclonal Antibody, Unconjugated Clone MEM-174
<b>Catalog Number</b>	853.163.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	IP, FC
<b>Immunogen</b>	Human granulocytes
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Concentration of 1mg/ml in phosphate buffered saline containing 15mM sodium azide, approx. pH 7.4
<b>Storage</b>	Store at +2-8°C. For long term storage freeze aliquots. Avoid freeze/thaw cycles.
<b>Note</b>	We recommend that each laboratory determine an optimum working titre for use in its particular application
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognizes CD11b antigen (MAC-1), a transmembrane protein type 1 mainly expressed on monocytes, granulocytes and NK-cells