

Product Datasheet

Creatine Kinase MB (CK-MB), 98% purity 190-24A

Product Name	Creatine Kinase MB (CK-MB), 98% purity
Catalog Number	190-24A
Alternative Names	Creatine Kinase MB, CKMB isoenzyme, CK MB, CKB
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.5 mg/mL (Siemens Centaur CP)
Storage	-20 °C
Purity	≥ 98% (electrophoresis)
CAS Number	9001-15-4
Source	Human Heart
Shelf Life	3 years
Biological Activity	Report U/mL (CK-NAC assay)
Specific Activity	Reported U/mg protein
Unit Definition	One unit will transfer one micromole of phosphate from Creatine Phosphate to ADP per minute at 37°C.