

## Product Datasheet

### Creatine Kinase MB (rCK-MB), Recombinant Human 190-24-R1

<b>Product Name</b>	Creatine Kinase MB (rCK-MB), Recombinant Human
<b>Catalog Number</b>	190-24-R1
<b>Tested Applications</b>	ELISA, WB
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Clear, colorless solution
<b>Concentration</b>	≥ 1.0 mg/mL (Coomassie)
<b>Storage</b>	-70 °C or below
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Recombinant
<b>Shelf Life</b>	2 years
<b>Biological Activity</b>	≥ 400 U/mL @ 37°C
<b>Unit Definition</b>	One unit will catalyze the transphosphorylation of one micromole of phosphate from creatine phosphate to ADP per minute at 37°C.