

Product Datasheet

Creatine Kinase MM (CKMM) Goat anti-Human Polyclonal Antibody 190-46

| | |
|----------------------------|---|
| Product Name | Creatine Kinase MM (CKMM) Goat anti-Human Polyclonal Antibody |
| Catalog Number | 190-46 |
| Species/Host | Goat |
| Tested Applications | ELISA, Blotting |
| Immunogen | High purity human Creatine Kinase MM (CKMM) |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid |
| Concentration | ≥ 1.0 mg/mL ($E^{0.1\%}_{280nm} = 1.4$) |
| Storage | -20 °C |
| Purity | Affinity Purified |
| Source | Polyclonal Antibody |
| Shelf Life | 2 years |