

# Product Datasheet

## Anti-Human CD15 Monoclonal Antibody, Unconjugated Clone MEM-158 853.183.020

Product Name	Anti-Human CD15 Monoclonal Antibody, Unconjugated Clone MEM-158
Catalog Number	853.183.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP
Immunogen	Human granulocytes
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. For longer storage freeze aliquots.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgM
Clonality	Monoclonal
Specificity	Reacts with CD15, a cell membrane molecule 3-fucosyl-N-acetyllactosamine (3-FAL) expressed on granulocytes, monocytes, macrophages, mast cells.