

Product Datasheet

Anti-Human CD22 Monoclonal Antibody, PE Conjugated Clone MEM-01 853.192.010

Product Name	Anti-Human CD22 Monoclonal Antibody, PE Conjugated Clone MEM-01
Catalog Number	853.192.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	Burkitt's lymphoma cell line Raji
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) containing 15 mM sodium azide and 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA)
Storage	Store at +2-8°C protected from light. DO NOT FREEZE.
Note	For flow cytometry we recommend using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD22 (BL-CAM), a 130 kDa type I transmembrane glycoprotein (immunoglobulin superfamily) expressed in the cytoplasm of pro-B and pre-B lymphocytes