

Product Datasheet

D-Dimer, anti-Human Monoclonal Antibody 200-14

Product Name D-Dimer, anti-Human Monoclonal Antibody

Catalog Number 200-14

Description Purified Monoclonal Antibody

Species/Host Mouse

Immunogen Purified, native human D-Dimer

Alternative Names anti-D-dimer

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration Reported in mg/mL

Storage -20 °C

Isotype IgG

Clonality Monoclonal

Purity Highly Purified

Shelf Life 2 years