Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Anti-Human CD22 Monoclonal Antibody, Unconjugated Clone MEM-01 853.193.020

Product Anti-Human CD22 Monoclonal Antibody, Unconjugated Clone

Name MEM-01

Catalog 853.193.020

Number

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested

WB, FC, IP

Applications Immunogen

Burkitt's lymphoma cell line Raji

Brand Diaclone

Form/AppearaPhosphate buffered saline (PBS) with 15 mM sodium azide,

approx. pH 7.4

Storage Store at +2-8°C. DO NOT FREEZE.

Note We recommend that each laboratory determine an optimum

working titre for use in its particular application

Isotype lgG1

Clonality Monoclonal

Specificity Reacts with CD22 (BL-CAM), a 130 kDa type I transmembrane

glycoprotein (immunoglobulin superfamily) expressed in the

cytoplasm of pro-B and pre-B lymphocytes